Installation manual

MANUEL D'INSTALLATION
INSTALLATIEHANDLEIDING
INSTALLATIONSANLEITUNG
MANUALE DI INSTALLAZIONE
MANUAL DE INSTALACIÓN
INSTRUKCJA INSTALACJI
PŘÍRUČKU PRO INSTALACI

Forest 3+

Contents

EN

Preface	4
1. This product	6
2. Scope of delivery	7
3. Assembly and adjustments	8
3.1. Tools	8
3.2. Adjusting the wheelchair to the user	9
3.3. Changing the battery	18
3.4. Changing wheels or tyres	19

Preface

To support you on the installation and repairs of this powered wheelchair, we provide you with this installation manual. Please read it carefully. If you still have questions after reading this manual, do not hesitate to contact Vermeiren.

The information in this manual applies to the following wheelchair(s):

- Forest 3+

Important note

Pictures of the product are used to clarify the instructions. Details of the depicted product may deviate from your product.

Information available

On our website <http://www.vermeiren.com/> you will always find the most recent version of the information in this manual. Please consult this website regularly for possible updates.

Visually impaired people can download the electronic version of this manual and have it read out by means of a text-to-speech software application.

	User manual For user and specialist dealer
	User manuals of applied operating console and battery charger For user and specialist dealer
	Installation instructions For specialist dealer
	Service manual for wheelchairs For specialist dealer

	Drawings of (spare) parts For specialist dealer
	EC declaration of conformity

1. This product

1. Backrest
2. Armrests
3. Arm pads
4. Operating console
5. Seat
6. Safety belt
7. Leg rests or footrests
8. Footplates

9. Steering wheels (front wheels)
10. Driving wheels (rear wheels)
11. Driving motors
12. Battery housing
13. Front light (optional)
14. Rear light (optional)
15. Anti-tipping
16. Identification plate

2. Scope of delivery

The following components are part of the delivery:

- Footrests
- Batteries
- Battery charger
- Tools
- User manual of:
 - wheelchair
 - Operating console
 - Battery charger
- Options

Before use, check if everything is included and that no products are damaged.

Please note that the basic configuration may differ in different European countries. Contact your specialist dealer.

3. Assembly and adjustments

⚠ CAUTION

Risk of injuries or damage

- Only use the settings described in this manual.
 - Only use parts described in this manual.
 - Only use genuine spare or replacement parts approved by Vermeiren.
 - Prevent bystanders from accessing the working area.
 - Keep the (dis)assembling area clean and tidy, to prevent falling.
 - Collect all small parts in a box or bag.
 - While (dis)assembling and operating the wheelchair, make sure no objects or body parts get caught between moving parts. In particular, pay attention to the position of the cables.
 - Turn the wheelchair off before you perform any mechanical adjustments and make sure that the free-wheel lever is in braking position.
 - Make sure all screws are firmly secured before driving.
 - Software changes are only allowed if executed by skilled personnel, taking into account the user and environment (stability, acceleration, indoor/outdoor, etc.).
- i** Read the user manual of this wheelchair. Be aware of the technical details and limits of intended use. For more information, contact Vermeiren.

3.1. Tools

To (dis)assemble and adjust the wheelchair, the following tools are needed:

- Wrench no. 10 to no. 17
- Allen key n° 4
- Screwdriver

3.2. Adjusting the wheelchair to the user

3.2.1. Seat inclination (manual adjustment of SE30)


CAUTION

Risk of injuries or damage

- The frame with seat is heavy to handle. Ask help of an assistant.
- Ask assistance to prevent the mounting frame with seat from falling during adjustment.
- Do not use the lowest seat inclination combined with 260mm front wheels and B06 footrests. These settings can lead to collision between the front wheels and footrests.

The seat angle can be manually adjusted according to the following steps:


1. Fold the footplate up or if possible remove the footrests.
2. Remove the screw (A) that connects the seat frame with the lower frame. Fold the seat frame backwards.
3. Check that the seat frame stays backwards or ask someone to hold the seat frame to the back.
4. Remove the screw (B).
5. Move the inner square tube (C) until the desired seat angle is reached (6 positions: 2° - 6° - 9,5° - 13° - 17° - 21°).
6. Retighten the screw (B).
7. Fold the seat frame back to the front.
8. Retighten the screw (A) to attachment (D).
9. Turn the footplate in place or install the footrests.

3.2.2. Seat depth

The wheelchair is adjustable in 6 different seat depths over a range of 100 mm (steps of 20 mm).


1. Loosen the two M6 socket screws (25) on either side of the frame.
2. Slide seat frame (26) backwards or forwards to get the desired position. Make sure that inner and outer tube holes are in line.
3. Insert the 4 socket screws (25) and tighten the nuts firmly.

3.2.3. Backrest inclination (manual adjustment of SE34)

1. Remove the M6 socket screw (10).

- Slide the inner square tube (5) in/out to the desired backrest angle (5 positions: 2° - 9° - 16° - 24° - 32°).
- Retighten the socket screw (10).
- Check that the backrest frame is firmly secured.

3.2.4. Width between armrests

⚠ CAUTION

Risk of injuries

- The maximum width is indicated by two triangles on each side of the seat, see picture. The outer triangles may not be covered by the frame tube.

To adjust the width of the armrests:


- Loosen the two M8 socket screws (32) at the bottom of the seat frame, behind the securing hooks.
- Adjust the armrest to a comfortable width.
- Tighten the socket screws (32) firmly.
- Repeat this at the other side of the wheelchair.
- Make sure the armrests are positioned identically on both sides of the seat.

3.2.5. Armrest length

The telescopic armrest can be adjusted in 6 positions as follows:

1. If necessary, cut the cable ties that hold the power cable of the console.
2. Loosen the socket screw (33).
3. Slide the inner tube (34) in/out to the desired length.
4. Insert socket screw (33) and tighten properly.
5. Repeat this for the second armrest.
6. Make sure both armrests have the same length.
7. Make sure that both armrests are firmly secured.
8. Tighten the cable with cable ties.

3.2.6. Armrest lowest and highest angular position

⚠ CAUTION

Risk of injuries

- Always use two socket bolts and two adjustment holes to tighten the armrest. Never exceed the maximum angular adjustment (holes).

To change the end position of the turning angle:

1. Hold the armrest and remove the two socket screws (35).
2. Turn the armrest support at the toothed clamp (36).
3. Make sure that the holes of inner and outer clamp part match.
4. Insert the two socket bolts again and tighten.
5. Repeat this for the second armrest.
6. Make sure both armrests have the same end positions.
7. Make sure that both toothed clamps (36) are firmly secured.


3.2.7. Operating console

3.2.7.1. Mounting position left or right

1. Cut the cable ties and unplug the power cable from the operating console.
2. Hold the control with its support (3) while removing M6 socket screw (41).
3. Loosen hex nut (42) and demount shaft (43) from the armrest.
4. Mount shaft (43) in the other armrest and lock it with hex nut (42).
5. Mount the control support (3) with socket screw (41) on the shaft (43).
6. Make sure all connections are well secured.

7. Guide the cable and plug the cable in the operating console.
8. Tie the cable to the back of the wheelchair subframe and on the armrest with cable ties. Make sure that the cable is long enough to make the adjustment of armrest and console possible.

3.2.7.2. Operating console height

The height of the operating console can be adjusted in 2+4 positions, or 6+4 positions depending on the type of shaft:


1. Hold the control while removing M6 socket screw (41).
2. Move the control support (3) up/down on shaft (43) and mount the socket screw (41).
3. If necessary adjust the height of the console on shaft (45) in the same way and lock it with socket screw (44).
4. Make sure all connections are well secured.
5. Guide the cable and plug the cable in the operating console.

3.2.8. Footrest length B06

CAUTION


Risk of damage

- Make sure that the footrests won't touch the ground. Keep a minimum safe distance of 60 mm above the ground.

1. Remove the screw (39).
2. Adjust the footrest to a comfortable length.
3. Retighten the socket screw firmly.
4. Repeat this for the second footrest. Make sure the footrests are positioned identically.

3.2.9. Footrest width

1. Loosen the two screws (A) on the front of the seat frame.
2. Move the footrest in an horizontal direction until the correct seat with, by moving the square tubes (B) and (C) over each other (3 different positions: steps 15 mm on each side).
3. Retighten the screws (A) properly.
4. Repeat instructions above for the other footrest.

3.2.10. Leg rest BZ7-BZ8


3.2.10.1. Leg rest length

⚠ CAUTION

Risk of injuries or damage


- The maximum length of the footrest is limited to the length indicated by a triangle on the front of the footrest tube, see picture. This maximum length adjustment is not possible in sitting position; the footplates will touch the ground.
- Make sure that the distance to the ground in sitting position is minimal

60 mm to avoid touching small obstacles and to prevent damage to the footrest.

To adjust the length of the leg rests:


1. Remove the socket screw (24) on the back of the footrest.
2. Adjust the leg rest to a comfortable length.
3. Retighten the socket screw firmly.
4. Repeat for the second leg rest. Make sure that the leg rests are positioned identically.

3.2.10.2. Footplate angle


To adjust the angle of the footplate at its toothed clamp:

1. Loosen the socket screw (25) on the leg rest (26).
2. Turn the footplate to the desired angle (80° - 100°). The dashes (30) on the toothed clamp indicate the angular positions.
3. Retighten the socket screw (25).
4. Repeat for the second leg rest.


3.2.10.3.Calfrest height

1. Turn the calfrest outwards.
2. Loosen the socket screw (30) that clamps the calfrest support on the footrest.
3. Slide the calfrest up/down to a comfortable position (step less).
4. Retighten the socket screw (30) with nut.
5. Repeat this for the second calf rest.
6. Make sure that the calf rests are positioned identically.


3.2.10.4.Calfrest depth

1. Loosen the socket screw (31) on the calfrest support (32).
2. Adjust the depth of the calfrest to a comfortable position (7 positions available).
3. Retighten the socket screw.
4. Repeat this for the second calf rest.
5. Make sure that the calf rests are positioned identically.

3.2.10.5.Width between calfrests

1. Loosen the two socket screws (33) on the back of the calf rest.
2. Adjust the width to a comfortable position (2 positions available). Do not use the hole (34) in the middle!
3. Retighten the two socket screws (33).
4. Make sure that the calf rests are positioned symmetrically.

3.3. Changing the battery


⚠ CAUTION

Risk of injuries or damage

- The batteries must only be replaced by trained personnel.
- Do not come into contact with the acid from the batteries.
- Replace the batteries in a well-ventilated room.
- Make sure that the ventilation holes of the battery housing are free from any contaminations.
- First read the charging instructions in the user manual before charging the batteries.
- Do not place fingers between the batteries and the battery housing.

i Replace both batteries at the same time and with the same type as the ones used in the wheelchair. If this is not possible, contact the nearest Vermeiren facility.

1. Switch OFF the wheelchair.
2. Take off the front battery cover by removing the two knobs on the top and the knob on the bottom.
3. Move aside the plate with power module by loosening and removing the screws in the top corners. You can now reach the batteries.
4. Switch ON the wheelchair to check if it works properly.

3.4. Changing wheels or tyres

CAUTION

Risk of injuries or damage


- Make sure that all air is released from the tyre, before removing it.
 - Use a wooden block and jack to lift the wheelchair. Always place the jack under the lower frame, never under any other parts.
 - Make sure that no body parts, the inner tube nor any other objects are pinched between tyre and rim.
 - Handle carefully to prevent damage to the rim.
1. Replace the inner tube, and if necessary also the tyre.
 2. Inflate the tyre using inflating equipment which complies to regulations and with a pressure indication in “bar”.
 3. Make sure that the pressure is correct, see the pressure indication on the tyres.

3.4.1. Driving wheels

-  If a tyre needs to be replaced, first check the type of tyre that has been fitted. In case of a tyre with line profile, do not switch left and right tyre.

1. Lift the driving wheel from the ground using a jack.
2. Remove the socket screws (24) and take the wheel from the flange.
3. Open tyre valve and push it gently to release the air.
4. Loosen the 5 screws that connect both rim halves (on the inner rim half).
5. Separate the rim halves and remove the inner tube from the tyre.
6. Inflate the new or repaired inner tube partly.
7. Insert the inner tube into the tyre.
8. Place the two rim halves back on the tyre. Make sure that the tube valve pops out of the opening in the rim.
9. Connect the rim halves with the 5 socket screws. Tighten by hand.
10. Inflate the tyre at the correct pressure, see the indication on the tyre.
11. Remove grease and other contamination from the thread of the 5 flange socket screws. Apply screw adhesive (ex. Loctite) to these screws.
12. Mount the wheel on the flange.

13. Lower the jack.
14. Make sure the socket bolts are firmly secured.

3.4.2. Steering wheels

3.4.2.1. Dismantling

1. Loosen the screw on the steering wheel axle (25) and remove it from the steering wheel fork.
2. Let the air out of the steering wheel by lightly pressing the pressure pin on the valve.
3. Loosen the 5 screws that hold the split rim together (26).
4. Now separate the rim sides.

3.4.2.2. Assembly

1. Insert the partly-filled inner tube into the tyre.
2. Connect the two sides of the rim through the tyres and use the 5 connecting screws to screw the rim together (26).
3. Make sure that the valve juts out of the valve opening in the rim.
4. Put the wheel back into the front wheel fork (25) and inflate the wheel.

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All rights reserved, including translation.