

* - option
 ** - not depicted

EP00267 / P


EV ILLICO				Plan No.
				1/2 EP00267 / P
Scale: 1:24	Drafted	Verified	Confirmed	
Date	02/11/2017	04/09/2024	04/09/2024	
Initials	DB	PM	AZ	

Material			
No.	Plan/Article	Description	Quantity
1	EP00247	LOWER FRAME	1
2*	EP00357	CENTRAL BRAKE	1
	EP00359	CENTRAL BRAKE FOR ONE SIDE	1
3	EP00256	LEFT SIDE RAILS	1
	EP00280	LEFT SIDE RAILS 969	1
4	EP00257	RIGHT SIDE RAILS	1
	EP00280	RIGHT SIDE RAILS 969	1
5	EP00251	FOOT SIDE WITH ABS PLATFORM	1
	EP00276	FOOT SIDE WITH ABS WITHOUT MOTOR	1
	EP00259	FOOT SIDE WITH METAL LATHS	1
	EP00277	FOOT SIDE WITH METAL LATHS WITHOUT MOTOR	1
	EP00367	FOOT SIDE WITH METAL LATHS MOTOR 2ND POS.	1
	EP00368	FOOT SIDE WITH ABS MOTOR 2ND POSITION	1
	E00026603	FOOT SIDE WITH WOODEN LATHS	1
6	EP00260	HEAD SIDE WITH ABS PLATFORM	1
	EP00261	HEAD SIDE WITH METAL LATHS	1
	E00026487	HEAD SIDE WITH WOODEN LATHS	1
7	EP00262	HEAD PANEL	1
	EP00272	HEAD PANEL 965	1
	EP00274	UPHOLSTERY HEAD PANEL	1
	EP00358	HEAD PANEL FULL	1
	EP00519	HAED SIDE PANEL / ILLICO VINTERNO	1
8	EP00263	FOOT PANEL	1
	EP00272	FOOT PANEL 965	1
	EP00273	UPHOLSTERY FOOT PANEL	1
	EP00358	FOOT PANEL FULL	1
	EP00520	FOOT SIDE PANEL / ILLICO VINTERNO	1
9*	EP00366	BUMPER FOR STANDARD PANEL	1
10	EP00268	SELF RAISING AID	1
11	EP00269	TRANSPORT SYSTEM	1
12*	EP00386	HORIZONTAL KIT FOR TRANSPORT SYSTEM	1
13*	EP00371	HORIZONTAL TRANSPORT SYSTEM	1
14*	EP00279	BED ELONGATION	1
15*	EP00411	VERTICAL TRANSPORT SYSTEM	1
16*	EP00473	WIDTH EXTENSION KIT	1
17*	EP00493	SPLIT SIDE RAILS	1
18**	EP00278	STICKERS	1

Material			
No.	Plan/Article	Description	Quantity
19**	U000369	CARTON FOR BED WITH STANDARD PANEL	1
	U000367	CARTON FOR BED WITH PANEL 965	1
	U000386	CARTON FOR BED WITH CENTR.BR. STAND. PANEL	1
	U000387	CARTON FOR BED WITH CENTR.BR. PANEL 965	1
	U000424	CARTON FOR BED WITH STANDARD PANEL VINTERNO	1
	U000425	CARTON FOR BED WITH PANEL 965 VINTERNO	1
	U000426	CARTON FOR BED WITH CENTR.BR. STAND. PANEL VINT.	1
	U000427	CARTON FOR BED WITH CENTR.BR. PANEL 965 VINTERNO	1
20**	2993398	CARTON CORNER	4
21**	U000430	CARTON BOX FOR SPLIT SIDERAILS	1
		CARTON BOX FOR WIDTH EXTENSION	1

EP00267 / P

* - option
 ** - not depicted

EV ILLICO				Plan No.
				2/2 EP00267 / P
Scale: ---	Drafted	Verified	Confirmed	
Date	02/11/2017	04/09/2024	04/09/2024	
Initials	DB	PM	AZ	