



| Material | | | |
|----------|------------------|-------------|----------|
| No. | Plan | Article | Quantity |
| 1 | YP08149 | 2608454 (L) | 1 |
| | YP08148 | 2608354 (R) | 1 |
| 2 | E001038 | 1913465 (R) | 1 |
| | | 1913464 (L) | 1 |
| 3* | XP03265(short) | 2996861 | 1 |
| 4 | X005565 | 1912287 | 1 |
| 5 | M6x45 DIN912 | 1003109 | 2 |
| 6 | M6 DIN986(Black) | 2994059 | 2 |
| 7 | M5X20 DIN912 | 1915452 | 3 |
| 8 | M5 DIN985 | 2994006 | 4 |
| 9 | XP03080 | 2996656 | 1 |

* Option

Only right side depicted

EP00460 / D

| | | | | |
|-------------------------|------------|------------|------------|-------------|
| INOVYS 2/35 BI5 ARMREST | | | | Plan No. |
| | | | | EP00460 / D |
| Scale: 3:10 | Drafted | Verified | Confirmed | |
| Date | 24/06/2021 | 12/08/2024 | 12/08/2024 | |
| Initials | AS | PM | AS | |