



Use WEICONLOCK AN 302-43.

Use WEICONLOCK AN 302-50.

Detail B
Scale: 3:5

Drawing shows only right side

EP00373 / C

BACKREST ADJUSTABLE ANGLE AND ADJUSTABLE HIGHT TRIGO Page 1/2				Plan No. EP00373/ C
Scale: 3:10	Drafted	Verified	Confirmed	
Date	13/02/2019	20/02/2020	20/02/2020	
Initials	AS	FD	TK	

Material			
No.	Plan	Article	Quantity
1	YP07083(35-45)	2745665	1
	YP07084(40-50)	2745765	1
2	YP07064(25-30)	2744965	1
	YP07065 (30-35)	2745065	1
	YP07066 (35-40)	2745165	1
	YP07067 (40-45)	2745265	1
	YP07068 (45-50)	2745365	1
3	XP02061	GP14700	1
4	XP01862	GP13200	1
5	M4x12	2994126	2
6	XP01632	GP09399	1
7	XP01861	GP12900	1
8	M5x12 DIN 7991	2994076	2
9	M6x30	1992721	1
10	M6	2994059	1
11	WASHER_M6	2994068	1
* 12	Hand grip GR-KS	2993809	1
	GR-ERGO-S(120/80)	2993810	1
13	Y004353	F000411	1

* Optional (depending on the length)

EP00373 / C

BACKREST ADJUSTABLE ANGLE AND ADJUSTABLE HIGHT TRIGO Page 2/2				Plan No. EP00373/ C
Scale: 3:10	Drafted	Verified	Confirmed	
Date	13/02/2019	20/02/2020	20/02/2020	
Initials	AS	FD	TK	