



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

TRIGO 30-36 CROSS			
Scale:	1:3	Drafted	Verified
Date	24.09.2018	22.01.2026	22.01.2026
Initials	PM	PM	AS

Plan No. EP00337/I

Material			
No.	Plan	Article	Quantity
Width 30cm			
1	YP06856	2757165 (32cm)	1
	YP06855	2757065 (34cm)	1
	YP06854	2756965 (36cm)	1
2	YP06853	2756865 (32cm)	1
	YP06852	2756765 (34cm)	1
	YP06851	2756665 (36cm)	1
3	YP06124	2727565	2
Width 32cm			
1	YP06859	2758265 (32cm)	1
	YP06858	2758165 (34cm)	1
	YP06857	2758065 (36cm)	1
2	YP06862	2757965 (32cm)	1
	YP06861	2757865 (34cm)	1
	YP06860	2757765 (36cm)	1
3	YP06125	2727865	2
Width 34cm			
1	YP06865	2758565 (32cm)	1
	YP06864	2758465 (34cm)	1
	YP06863	2758365 (36cm)	1
2	YP06868	2758865 (32cm)	1
	YP06867	2758765 (34cm)	1
	YP06866	2758665 (36cm)	1
3	YP06125	2727865	2
Width 36cm			
1	YP06871	2759165 (32cm)	1
	YP06870	2759065 (34cm)	1
	YP06869	2758965 (36cm)	1
2	YP06874	2759465 (32cm)	1
	YP06873	2759365 (34cm)	1
	YP06872	2759265 (36cm)	1
3	YP05691	2719865	2

Material			
No.	Plan	Article	Quantity
4	YP06152 (raw)	27943	2
5	XP01851	2993570	1
6	XP01439	2993524	8
7	P00030394	2995366	4
8	X002136	1992414	4
9	CLOSED NUT M6	2994059	2
10	M8x80	2994121	1
11	Nut M8	1996213	1
12	M6 DIN125	2994068	2
13	Screw D3.9X9.5	2992349	4
14	M6x45 DIN912	2994058	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

TRIGO 30-36 CROSS			
Scale: 1:3	Drafted	Verified	Confirmed
Date	24.09.2018	22.01.2026	22.01.2026
Initials	PM	PM	AS

Plan No.
EP00337/I
