



| Material |             |         |     |
|----------|-------------|---------|-----|
| No.      | Plan        | Article | Qty |
| 1        | X005003     | 1907931 | 1   |
| 2        | 6002 2Z     | 3990423 | 1   |
| 3        | ---         | 3992846 | 1   |
| 3A       | XC00114     | 3990424 | 1   |
| 3B       | 6002 2Z     | 3990423 | 1   |
| 4        | ---         | 3992875 | 1   |
| 5        | M6x8 DIN916 | 2996905 | 2   |

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

3990495

PLASTIC FORK FOR D100 WITH BEARINGS

Plan No.

E00074488/A.1

|          |            |            |            |            |
|----------|------------|------------|------------|------------|
| Scale:   | 3:4        | Drafted    | Verified   | Confirmed  |
| Date     | 31/01/2024 | 31/01/2024 | 31/01/2024 | 31/01/2024 |
| Initials | SD         | SD         | PJ         |            |

