

Use WEICONLOCK AN 302-43

Use WEICONLOCK AN 302-43

Detail B
Scale: 3:5

Drawing shows only right side

BACKREST ADJUSTABLE ANGLE TRIGO

Page 1/2

Plan No.

EP00356/ G

Scale: 3:10

Drafted

Verified

Confirmed

Date

29/10/2018

10/09/2020

10/09/2020

Initials

PM


MK

TK



EP00356 / G

Material			
No.	Plan	Article	Quantity
* 1	YP06081(35)	2793265	1
	YP06082(37)	2793365	1
	YP06083(39)	2793465	1
	YP06084(41)	2793565	1
	YP06085(43)	2793665	1
	YP06086(45)	2793765	1
	YP06087(47)	2793865	1
	YP06088(49)	2793965	1
	YP06089(51)	2794065	1
2	XP01632	GP09399	1
3	XP02061	GP14700	1
4	XP01862	GP13200	1
5	M4x12	2994126	2
6	---	---	---
7	XP01861	GP12900	1
8	M5x12 DIN 7991	2994076	2
9	CAP ϕ 22	1990648	1
** 10	Fold. push handle	EP00246	1
** 11	M6x10	2994127	4
12	YP08333(35)	2613865	1
	YP08334(37)	2613965	1
	YP08335(39)	2614065	1
	YP08336(41)	2614165	1
	YP08337(43)	2614265	1
	YP08338(45)	2614365	1
	YP08339(47)	2614465	1
	YP08340(49)	2614565	1
	YP08341(51)	2614665	1
	13	Hand grip GR-KS	2993809
GR-ERGO-S		2993810	1



Plan	Lenght "X"
YP08333(35)	380mm
YP08334(37)	400mm
YP08335(39)	420mm
YP08336(41)	440mm
YP08337(43)	460mm
YP08338(45)	480mm
YP08339(47)	500mm
YP08340(49)	520mm
YP08341(51)	540mm

Plan	Lenght "X"
YP06081(35)	320mm
YP06082(37)	340mm
YP06083(39)	360mm
YP06084(41)	380mm
YP06085(43)	400mm
YP06086(45)	420mm
YP06087(47)	440mm
YP06088(49)	460mm
YP06089(51)	480mm

* Gdy uzywamy z raczka skladana(10,11) wybrac 2 rozmiary krotsza rure.

* When using a folded handle (10,11), select 2 sizes shorter tube.

** Option.

** Opcja.

EP00356 / G

BACKREST ADJUSTABLE ANGLE TRIGO				Plan No.
Page 2/2				EP00356/ G
Scale: 3:10	Drafted	Verified	Confirmed	
Date	29/10/2018	10/09/2020	10/09/2020	
Initials	PM	MK	TK	