



Material			
No.	Plan	Article	Quantity
1	YP08202 (R)	2610054	1
	YP08216 (L)	2610154	1
2	YP08206	2610300	2
3	YP08207	2610254	2
4	XP02978	2996607	2
5	XP02979	2996608	2
6	XP02093	GP15800	2
7	XP02466	GP21354	4
8	X004747	1906057	4
9	M5x25	2994906	2
10	M6x30	2994129	2
11	M6x55	2994376	4
12	M6x16	2994067	4
13	Nut M6	2994063	6
14	Nut M5	2994006	2
15	Washer M8	2994914	4
16	Tube cap D22	1990648	4

Only left side depicted, quantities for pair

EP00509 / B

ANTITIPPING TRIGO TYPE CH				Plan No.
				EP00509/ B
Scale: 1:3	Drafted	Verified	Confirmed	
Date	30/03/2020	23/10/2020	23/10/2020	
Initials	PM	AS	TK	