



no	part nu	part	QT	no	part nu	part	QT
1	1920656	MAIN FRAME		1	36		
2	1918862	FRONT FORK		1	37		
3					38	1920657 BATTERY CONNECTOR HOLDER	1
4	1920340	FRONT AXLE L			39	1920658 BATTERY CONNECTOR	1
5	1920341	FRONT AXLE R			40		
6					41		
7					42		
7					43		
7					44		

7				46	1918857	BALL JOINT LE	
8				45	1918858	BALL JOINT RIGHT	
8				47			
8				48			
8				49	1918859	CONNECTION ROD	
9			1	50			
10				51	1918851	FOOT PLATE	1
11				52	1920443	SET COLOR PANEL FRONT green	
				52	1920444	SET COLOR PANEL FRONT orange	
				52	1920445	SET COLOR PANEL FRONT Blue	
				52	1920446	SET COLOR PANEL FRONT pink	
			1	53			
13				54	1920634	WASHER M4 NF E25-513	
14				55			
15			2	56			
16				57	1918861	FRONT CONTACT POINT HOLDER	1
17				58			
18			2	59	1918860	CONTACT POINT FRONT	4
19				60	1920635	BOLT CAP	6
20				61			
21				62			
22	1918852	FRONT COVER	1				
23	1919090	POWER MODULE	1				
24	1918856	STEERING AXLE ONE	1				
25	1920633	BOLT M4X20 DIN 7985	13				
26							
27							
28							
29							
30							
31							
32	1920334	FRONTWHEEL ONE AIR	2				
33	1920632	FRONT BUMPER	1				
34							
35							




no	part nu	part	QT		part nu	part	QT
0	1920394	TIRE	4				
0	1001263	INNER TIRE	4				
1					17		
2					18		
3					19		
3	1918887	DIFF.+MOTOR+REM	1	20			
3	1920395	DIFF	1				
not depicted	1922014	BUSH TRANSMISION .	1				

3	1920396	MOTOR	1			
3	1920397	MOTORBRAKE	1			
4				21		
4				22		
4				23	1920342	REAR WHEEL LU 2
4				24		
5				25		
6				26		
7				27	1918889	REAR CONTACT POINT HOLDER 1
8				28		
9				30	1918890	CONTACT POINT REAR 4
10	1921454	REAR RED REFLECTOR	2			
11	1918882	REAR COVER	1	31		
12	1918886	HANDLE REAR	1	32		
13	1918853	COVER MIDDLE	1	33		
14	1920439	COLOR PANEL SET L/R green		34		
14	1920440	COLOR PANEL SET L/R orange				
14	1920441	COLOR PANEL SET L/R Pink				
14	1920442	COLOR PANEL SET L/R blue				
15				35		
16				36		




no	part nu	part	QT				
1							
2							
3	1921320	outer knob angle adjust	1				
4	1921319	middle knob angle adjust	1				
5	1921318	inner knob angle adjust	1				
6	1921322	bolt M5 X 16					
7	1921321	spring knob angle adjust	1				
8							
9							
10							
11	1921317	bolt angle adjust	1				
0		Charger	1				




no	part nu	part	QT				
0	1920391	BATTERYBOX 25 AHCOMPLETE	1				
0	1920392	BATTERY BOX 25AH WITHOUT BATTERIES	1				
1							
2							
3	1920531	MAIN FUSE					
4							
5							
6							
7							
8							

9						
10						
11						
12						
13						
14						
15	1920301	BATTERY 25 A/H		2		
16						
17						
18						
0	1919798	CHARGER		1		




no	part nu	part	QT			
0						
1						
2	1918898	SEAT TUBE	1			
3						
4	1918895	ARM PAD	2			
5	1920266	BOLT ARMREST FIXATION	2			
6	1918899	CHAIR	1			
6.1	1920640	BACK	1			
6.2	1920641	SEAT				
7	1920645	SEAT HANDLE	1			
8						
8						
9						
10						
11						
12						
13						
14						
15						
16						
17						
18	1918903	ARMREST LOWER LEFT	1			
19	1918907	ARMREST LOWER RIGHT	1			
20						
21	1918931	ARMRESR UPPER PART	2			
22						
23						
24						
	1920830	ARM REST COMPLETE L	1			
	1920831	ARM REST COMPLETE R	1			

1921438 SEAT PLATE COMPL W/A tube



no	part nu	part	QT
1			
2	1918911	SPEED DIAL ASSEMBLY	1
3	1918918	CHARGER SOCKET	1
4			
5	1918910	SPEED LEVER	1
6			
7	1918913	STEERING COVER BACK	1
8	1921463	FRONT LIGHT COVER	1
9	1918914	STEERING COVER FRONT	1
10	1919437	CONTACT IGNITION	1
11	1920423	HORN SWITCH	1

12	1918912	BATTERY INDICATOR	1
13	1918909	POT SPEED	1
14			
15			
16			
17			
18			
19			
20			
21	1920502	ASSEMBLY PLATE BASKET MIDDLE	1
22	1920510	ASSEMBLY PLATE BASKET REAR	1
23	1920501	ASSEMBLY PLATE BASKET FRONT	1
	1920832	BASKET HOLDER COMPLETE	
24	1920511	BOLT M5X25	2
25	1920503	WASHER	4
26	1920505	BOLT M5X20	4
27			
28	1918908	IGNITION KEY	1
29	1920509	BASKET	1
30			
31			
32	1920629	FRONT LIGHT	1
33			
34			
35			
	1918977	MIRROR LEFT	1
	1918978	MIRROR RICHT	1
37	1920628	HORN	1



no	part nu	part	QT
1	1918918	CHARGER SOCKET	1
2			
3			
4			
5			
6			
7	1918976	MAIN CABLE	1
8			
9	1918909	POT. SPEED	1
10	1921636	CABLE+ FUSE	1
11			
12			
13			
14			
15			
16			