



Material			
No.	Plan	Article	Quantity
1	YP05905(R,)	2726665	1
	YP05807(L,)	2725465	1
	YP07485(R)short	2770665	1
	YP07486(L)short	2770765	1
	YP06329(R)	2797965	1
2	YP06330(L)	2797865	1
	YP07490(R)short	2770965	1
	YP07491(L)short	2770865	1
	YP02372(R)S	1899865	1
	YP07780(R)M	2781965	1
3	YP08105(R)L	2607865	1
	YP02373(L)S	1899765	1
	YP07781(L)M	2782065	1
	YP08106(L)L	2607965	1
	XP00615	2990920	1
4	XP01843(Φ22)	2993563	1
	XP00697(Φ25.4)	2990919	1
5	XP00335	2990264	1
6	XP01560	2993561	1
7	---	2990351	0,4m
8	CAP_D22	1990648	1
9	M6X45	2994125	2
10	M6x40	2994123	1
11	M6x25	2994062	1
12	M6	2994063	3

EP00270 / D

TUBULAR ARMREST TRIGO				Plan No.
				EP00270 / D
Scale: 1:1	Drafted	Verified	Confirmed	
Date	28/06/2018	16/10/2019	16/10/2019	
Initials	AS	MK	TK	