



| Material | | | |
|----------|------------------------|---------|----------|
| No. | Plan | Article | Quantity |
| 1 | X002976 | 1995343 | 2 |
| 2A | P00065423 | 2684165 | 1 |
| 2B | P00071838 | 2684265 | 1 |
| 3 | Y004693 | 1518204 | 2 |
| 4 | Y004751 | 1521604 | 2 |
| 5 | M6x40 ISO7380 | 1006111 | 4 |
| 6 | M6 DIN7349 | 1003321 | 1 |
| 7 | M6x25 | 1903394 | 1 |
| 8 | BLACK WHEELS SET (L+R) | 1901208 | 1 |
| | GRAY WHEELS SET (L+R) | 1901745 | 1 |
| 9 | --- | --- | - |
| 10 | --- | 1994568 | 4 |
| 11 | DIN471 | 1995259 | 4 |
| 12A | P00064457 | 2683565 | 1 |
| 12B | P00071833 | 2683665 | 1 |
| 12C | P00071837 | 2684065 | 1 |
| 13A | X005892 | G103100 | 1 |
| 13B | X005891 | G103000 | 1 |
| 16 | M6 | 1008492 | 4 |
| 17 | XP03058 | GP23399 | 2 |
| 18 | P00064455 | GP29965 | 2 |
| 19 | D3x26 DIN 1481 | 1908821 | 2 |
| 20 | 5x5x10 DIN6885A | 2998153 | 2 |

Part 18

Protect socket/hole.
Free from painting
on both side.

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

| | | | | |
|-------------------------|------------|------------|------------|-----------|
| DRIVE SYSTEM HEM2 TRIGO | | | | Plan No. |
| | | | | EP00521/G |
| Scale: 1:6 | Drafted | Verified | Confirmed | |
| Date | 03.07.2020 | 26.01.2026 | 26.01.2026 | |
| Initials | AS | PM | AS | |

2B-12C- 13A Seat widths 38-40
 2A-12A - 13A Seat widths 42-44
 2A-12B - 13B Seat widths 46-48-50