


NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

VERSO B74 DRUM BRAKE WHEELS			
Scale:	1:4	Drafted	Verified
Date	31/07/2025	31/07/2025	31/07/2025
Initials	SD	SD	TS

Plan No.
E00113569/A1.1





Material			
No.	Plan	Article	Qty
1	P00027229	2643800	2
2	P00099263	1922329	2
3	P00099264	1922330	2
4	P00100842	1922331	2
6	X003986	1902199	2
7	X005514	1906244	2
8	Left Set	3990705	1
	Right set	3990708	1
9	X005545	2019900	2
10	Y017346	17480..	2
11	D25.4	3990707	2
12	D3x16 DIN1481	1003791	2
13	M6x16 ISO7380	1008459	2
14	24''' Black PU	1914616	2
15	M6x50 DIN912	2995258	2
16	M6 DIN986	2994059	2
17	M6 DIN9021	2994908	2
18	E00023606	1918195 (L)	1
	E00023606	1918196 (R)	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

VERSO B74 DRUM BRAKE WHEELS

Plan No.

E00113569/A1.1

Scale:	1:2	Drafted	Verified	Confirmed
Date	31/07/2025	31/07/2025	31/07/2025	31/07/2025
Initials	SD	SD	TS	

