



Detail A
Scale: 1:3

E000712 / F

NAVIX SEAT FRAME L35 ELECTRIC				Plan No.
Sheet 1/2				E000712 / F
Scale: 1:6	Drafted	Verified	Confirmed	
Date	16/08/2011	26/11/2015	26/11/2015	
Initials	JG	SD	PJ	

Material							
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	X001016	1001682	2	12	M8x12 DIN916	1004885	6
2	X005019	1909238	2	13	M6x12 ISO7380	1007000	2
3	X004804	1754800	2	14	---	1906790	2
4	---	1907467 (40)	1	15	---	1906791	2
	---	1907468 (45)	1	16	---	1906792	2
	---	1907469 (50)	1	17	M10 DIN985	1001683	1
5	Y017305	1751954 (40)	1	18	M8x16 ISO7380	1991858	4
	Y017285	1748754 (45)	1	19	---	1992835	3
	Y017282	1751354 (50)	1	20	M4x8 DIN7984	1996530	1
6	Y016954	1752354 (40)	1	21	M8x25 ISO7380	1990501	2
	Y017347	1750154 (45)	1	22	M10x35 ISO7380	1009442	1
	Y016960	1752454 (50)	1	23	M8x40 ISO7380	1008739	3
7	Y017783	1755454 (40)	1	24	X003111	1995853	2
	Y017784	1755554 (45)	1	25	M8	1902024	2
	Y017785	1755654 (50)	1	26	X005383	1911710	1
8	Y017331	1748954	1	27	X005384	2014400	1
9	Y017643	1754754	1	28	Y021834	2017254	1
10	---	1992805	4	29	M8x30 ISO7380	1992717	1
11	---	1001612	4	30	M8x35 DIN933	1001444	1
				31	M8 DIN934	1001671	1

E000712 / F

NAVIX SEAT FRAME L35 ELECTRIC				Plan No.
Sheet 2/2				E000712 / F
Scale: 1:6	Drafted	Verified	Confirmed	
Date	16/08/2011	26/11/2015	26/11/2015	
Initials	JG	SD	PJ	