

Material			
No.	Plan	Article	Qty
1	P00068432	2086754	2
2	P00085567	2092532	2
3	---	3990706	2
3A	D22	1994012	2
3B	---	1993218	2x0,75m
3C	---	1004094	2x0,95m
4	X000185	G000790	2
5	X002356	1502800	2
6	Left side	3990705	1
	Right side	3990708	1
6A	X002437	1992699	2
6B	M16 DIN17222	1911653	2
6C	X005816	1914818	2
6D	X005815	1914817	2
6E	---	1906205	2
9	XC00017	1C11600	2
10A	Y024006	F002907	2
10B	---	3990041	2
12	M6x50 DIN7991	1001486	2
13	M6x45 DIN912	1003109	2
15	M6x20 ISO7380	1006999	2
16	M6 DIN986	1008492	4
17	24"	1914616	2
18*	EC00510	3998373 (L)	1
	EC00511	3998374 (R)	1

* Brakes not depicted
10A for use with B90
10B for use without B90

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

D200N B74 DRUM BRAKE WHEELS				Plan No.
				E00087898/B.1
Scale:	1:4	Drafted	Verified	Confirmed
Date	30/08/2024	10/10/2024	10/10/2024	
Initials	SD	DVB	PJ	

